GIMA

**ASPIRATORE SUPER VEGA
ASPIRATOR SUPER VEGA
ASPIRATEUR SUPER VEGA
ABSAUGER SUPER VEGA
ASPIRADOR SUPER VEGA**

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SUPER VEGA it's a device working 230V ~ / 50 Hz network electricity, to be used for the nasal aspiration, oral aspiration, tracheal aspiration of the body liquids (mucus, catarrh or blood) in the adult or in the children. Easily portable equipment designed for not continuous use 120 min. ON / 60 min. OFF. Thanks to this characteristics and to the rating that it has, this product is particularly suitable for hospital use, on the tracheotomized patients, minor surgical applications and post-operative therapy at home. Made of highly heat resistant, electrically insulated plastic material in conformity with the latest European safety standard, the product is supplied with a complete polycarbonate autoclavable jar with overflow valve and it is equipped with aspiration regulator and vacuum indicator located on the front panel.

GENERAL WARNINGS

READ INSTRUCTION MANUAL CAREFULLY BEFORE USE

ONLY HIGHLY QUALIFIED STAFF USE RESERVED

**THE INSTRUMENT MUST NOT DISASSEMBLED
FOR TECHNICAL SERVICE ALWAYS CONTACT GIMA**




IMPORTANT SAFETY RULES

1. On opening the packaging, check the integrity of the appliance, paying particular attention to the presence of damage to the plastic parts, which may make access possible to internal live parts and also to breakage and / or peeling of the power supply cable. **In these cases don't connect the plug to the electric socket. Carry out these controls before each use;**
2. before connecting the appliance always check that the electric data indicated on the data label and the type of plug used, correspond to those of the mains electricity to which it's to be connected;
3. If the plug supplied with the appliance is incompatible with the mains electricity socket, contact qualified staff for replacement of the plug with a suitable type. The use of simple or multiple and / or extension adapters is not generally recommended. Whenever their use is indispensable, use those in compliance with safety regulations, however paying attention not to exceed the maximum power supply limits, which are indicated on the adapters and extensions;
4. Respect the safety regulations indicated for electrical appliances and particularly:
 - Only use original accessories and components;
 - The device can be used only with the bacteriological filter;
 - Never immerge the appliance into water;
 - Position the appliance on flat stable surfaces;
 - Position the device in a way that the air inlets on the back aren't obstructed;
 - Never use the device in environments which have anaesthetic mixtures inflammable with air, oxygen or nitric oxide;
 - Don't touch the device with wet hands and always prevent the appliance coming into contact with liquids;
 - Keep off the reach of children or not capable people without supervision;
 - Don't leave the appliance connected to the power supply socket when not in use;
 - Don't pull the power supply cable to disconnect the plug remove the plug from the mains socket correctly;
 - Preserve and use the medical device in environments protected from atmospheric factors and at a distance from heat sources;
5. For repairs, exclusively contact GIMA technical service and request the use of original spare parts. Failure to comply with the above can jeopardise the safety of the device;
6. **This medical device must be destined exclusively for the use for which it has been designed and described in this manual.** Any different use must be considered incorrect and therefore dangerous; the manufacturer cannot be considered liable for damage caused by improper, incorrect and / or unreasonable use or if the appliance is used in electrical plants that are not in compliance with the regulations in force;
7. Particular precautions must be made concerning electromagnetic compatibility. The medical device must be installed and used according to information supplied with the accompanying documents;
8. Instrument and accessories discharging must be done following current law regulations in every country of use.

TECHNICAL CHARACTERISTICS

| | | |
|--|---------------------------|---------------------------------|
| Model | SUPER VEGA | |
| Typology (MDD 93/42/EEC) | Medical Device Class II a | |
| Classification UNI EN ISO 10079-1 | HIGH VACUUM / HIGH FLOW | |
| Main Voltage | 230 V ~ / 50 Hz | 110 V ~ / 60 Hz (no CE 0123) |
| Power consumption | 107 VA | 165 VA |
| Fuse | F 1 x 1.6 A 250 V | F 1 x 4 A 250 V |
| Maximum suction aspiration (without jar) | - 80kPa (- 0.80 bar) | |
| Minimum suction aspiration | Less -40kPa (- 0.40 bar) | |
| Maximum flow (without jar) | 40 l/min | |
| Weight | 3.6 Kg | |
| Dimension | 350 x 210 x 180 mm | |
| Functioning (to 35°C and 110 % operating voltage) | 120 min ON / 60 min OFF | |
| Working condition | Room temperature: | 5 ÷ 35 °C |
| | Room humidity percentage: | 30 ÷ 75 % RH |
| | Altitude: | 0 ÷ 2000m s.l.m |
| Conservation condition | Room temperature: | - 40 ÷ 70 °C |
| | Room humidity percentage: | 10 ÷ 100% RH |

SYMBOLS

| | |
|---|--|
|  | Class 2 isolation equipment |
|  | CE marking in conformity with EC directive 93/42/CEE Manufacturer : CA-MI snc Via Ugo La Malfa n°31 43010 Pilastro (PR) Italy |
|  | Warning, consult the instruction manual |
|  | Type B equipment |
|  | Fuse |
|  | To preserve in place coolness and dry land |
|  | Conservation temperature: - 40 ÷ 70° C |
|  | Alternate current |
|  | Mains frequency |
|  | ON |
|  | OFF |

| Guidance and manufacturer's declaration – electromagnetic Emissions (as request by regulation EN 60601-1-2:2001) | | |
|--|-------------------|---|
| The surgical aspirator SUPER VEGA is intended for use in the electromagnetic environment specified below. The customer or the user of the surgical aspirator SUPER VEGA should assure that it's used in such an environment. | | |
| Emissions test | Compliance | Electromagnetic environment - guidance |
| Irradiated / Conducted emissions CISPR11 | Group 1 | The surgical aspirator SUPER VEGA only used RF energy only for its internal functioning. Therefore, its RF emissions are very low and are not cause interference in proximity of any Electronic appliances. |
| Irradiated / Conducted emissions CISPR11 | Class [B] | The surgical aspirator SUPER VEGA can be used in all environments, including domestic and those connected directly to the public mains distribution that supplies power to environments used for domestic scopes. |
| Harmonic emissions IEC/EN 61000-3-2 | Class [A] | |
| Voltage fluctuations / flicker emissions IEC/EN 61000-3-3 | Complies | |

| Guidance and manufacturer's declaration – electromagnetic Immunity (as request by regulation EN 60601-1-2:2001) | | |
|--|---|---|
| The surgical aspirator SUPER VEGA is intended for use in the electromagnetic environment specified below. The customer or the user of the surgical aspirator SUPER VEGA should assure that it's used in such an environment. | | |
| Immunity Test | Compliance | Electromagnetic environment - guidance |
| Electrostatic discharge (ESD) IEC/EN 61000-4-2 | ± 6kV on contact ± 8kV in air | Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%. |
| Electrical fast transient / burst IEC/EN 61000-4-4 | ± 2kV power supply | Mains power quality should be that of a typical commercial environment or hospital. |
| Surge IEC/EN 61000-4-5 | ± 1kV differential mode | Mains power quality should be that of a typical commercial environment or hospital. |
| Loss of voltage, brief voltage interruptions and variations IEC/EN 61000-4-11 | 5%U _T for 0.5 cycle 40%U _T for 05 cycle 70%U _T for 25 cycle <5%U _T for 5 sec | Mains power quality should be that of a typical commercial environment or hospital. If the user of the surgical aspirator SUPER VEGA request that the appliance operates continuously, the use of a continuity unit is recommended. |
| Magnetic field IEC/EN 61000-4-8 | 3A/m | The power frequency magnetic field should be measured in the intended installation location to assure that it's sufficiently low. |
| Conducted Immunity IEC/EN 61000-4-6 | 3Vrms 150kHz to 80MHz (for appliances that aren't life - supporting) | - |
| Irradiated Conducted IEC/EN 61000-4-3 | 3V/m 80MHz to 2.5 GHz (for appliances that aren't life – equipment) | - |
| Nota U _T is the value of the power supply voltage | | |

ACCESSORIES SUPPLIES

| DESCRIPTION |
|--|
| COMPLETE ASPIRATION JAR |
| CONICAL FITTING |
| TUBES SET 6mm x 10mm (TRASPARENT SILICON) |
| ASPIRATION PROBE CH 20 |
| ANTIBACTERIAL FILTER |

The filter is produced with (PTFE) hydrophobic material witch prevents fluids entering the pneumatic circuit. When the filter is wet, it's not possible to use the unit therefore the filter should be changed immediately. In case of possible contamination or discolouration, change the filter immediately. Don't use the suction unit without the protection filter fitted. If the suction unit is used in an emergency or in a patient where the risk of contamination is not know the filter must be changed after each use. The suction probe must be changed at the end of each application. Available under request with different versions with two complete jar 2200cc.

CLEANING OF ACCESSORIES

Wearing disposable latex gloves, proceed to clean the accessories as follows:

1. Grip the jar with your hand and turn the cap in a counter-clockwise direction.
2. Remove the red cage covering the float located in the jar cap (conical part with 3 openings)
3. Once the cage is removed, the float with the red seal comes out of the jar cap.
4. **To disinfect the accessories, use either denatured alcohol or hypochlorate-based solutions, easily purchasable from chemist shops.**
5. Once cleaning has been completed, insert the float in the cage, with the seal facing towards the cage opening.
6. Insert the opening of the cage into the seat in the jar cap.
7. To facilitate the cleaning of the jar, fill it 1/3 full of ordinary water. The aspirated liquid, thus diluted, will then be easier to remove.

It is possible autoclave the accessories cover and jar: insert the parts in autoclave and to effect a cycle of sterilization with vapour to the temperature of 121°C (pressure relative 1 bar) having care to position turned upside-down the graduated jar (with turned fund toward the tall one).

After the sterilization and the cooling to the temperature the components' environment, to verify that this last don't result damaged; disassembled therefore the container for inhaled liquids following the inverse operations to the dismantlement.

The aspiration tubes can be sterilized on autoclave using a sterilization cycle at 120°C. The conical connector can be sterilized on autoclave using a sterilization cycle at 121°C. The device is ready for a new employment now.

During the sterilization operations the operator must use disposable gloves in latex, shirts and oportune mask not to enter contact with possible substances contaminating

DO NOT WASH, STERILIZE OR PUT IN AUTOCLAVE THE ANTIBACTERIAL FILTER

PERIODICAL MAINTENANCE CHECKS

The SUPER VEGA suction equipment does not need maintenance or lubrication. It is necessary to check functioning and instrument before every use. Unpack the instrument and **always check** integrity of plastic parts and feeding cable, they might have been damaged during previous use. Connect cable to electrical network and turn switch on. Close the aspirator outlet with your finger and with suction regulator in maximum vacuum position check that the vacuum indicators reaches -80kPa (-0.80 bar) minimum. Rotate the knob from right to left and check the aspiration regulating control. The vacuum indicator should go down -40kPa (-0.40 bar). Verify that loud noises are not present, these can indicate wrong functioning. Two protection fuses (**F 1 x 1.6 A 250V**) reachable from exterior and situated in the plug protects the instrument. For fuses replacing, always the type and the range.

| Type of fault | Cause | Remedy |
|---|---|--|
| 1. No aspiration | Jar Cap badly screwed down | Unscrewed the cap, then rescrew it correctly |
| 2. No aspiration | Lid seal not in its seat | Unscrew the cap and insert the seal properly in its seat |
| 3. The float doesn't close | If the cap has been washed, ensure that the float is not partially detached | Insert the float into it's place |
| 4. The float doesn't close | The float it's covered by dirty material | Unscrewed the cap, leave the and put in on autoclave |
| 5. Low suction | Foam inside the jar | Fill the jar to 1/3 full of ordinary water |
| 6. No aspiration due to flow leakage of mucus | Filter blocked | Replace filter |
| Faults 1 - 2 - 3 - 4 - 5 - 6 | None of the remedies has achieved the desired results | Contact the seller or GIMA After-sales Assistance Service |

If the overflow security system it's activated, don't proceed with the liquid aspiration. If the overflow security system doesn't work there are two cases:
 1° case – If the overflow security system doesn't work the aspiration will be stopped by the bacteriological filter who avoid the liquid penetration inside the device.
 2° case – If both the security system doesn't work, there is the possibility that liquid comes inside the device, in this case return the device to GIMA technical service.

BEFORE EVERY CHECKING OPERATION, IN CASE OF ANOMALIES OR BAD FUNCTIONING, PLEASE CONTACT GIMA TECHNICAL SERVICE. GIMA DOES NOT GIVE GUARANTEE IF INSTRUMENT, AFTER THE TECHNICAL SERVICE CHECKING, APPEARS TO BE TAMPERED.

IMPORTANT INFORMATION FOR CORRECT DISPOSAL OF THE PRODUCT IN ACCORDANCE WITH EC DIRECTIVE 2002/96/EC:

At the end of its working life, the product must not be disposed of as urban waste. It must be taken a special local authority differentiated waste collection or to a dealer providing this service. Disposing of a household appliance separately avoids possible negative consequences for the environment and health deriving from inappropriate disposal and enables the constituent materials to be recovered to obtain significant saving in energy and resources. As a reminder of the need to dispose of separately the product is marked with a crossed-out wheeled dustbin.

INSTRUCTION FOR USE

- Connect the short silicon tube **1** ,with filter antibacterial **8**, to the suction connector **2**. The other tube, with one end connected to the filter must be connected with the other end to jar's lid connector **3** where has been fixed the red float. When the 90% of the volume of the jar is reached there is the activation of the security float (the float close the aspiration connector on the jar) to avoid liquid penetration inside the device. The device must be used on a plan of horizontal operation.
- Connect the long silicon tube **4** to the other jar's lid connector **5**
- Connect the other end of the long silicon tube **4** to the probe plastic connector **6** then connect the suction probe to it.
- Connect the power cord to the device then connect the plug to the electrical mains supply.
- Push switch **7** on position **1** to start suction
- Once finished push switch on **0** position and unplug.
- Unscrew the jar's lid and fill the jar 1/3 full or ordinary water (this for an easy cleaning operations and an rapid reaching of the functionally vacuum) then rescrew the lid on the jar correctly.
- To extract the accessories and start with cleaning.

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